

Hydraulic to Electric Conversion for Midmark Tables

Stop Fixing Old Hydraulics!

We install New Electric Actuators in Midmark Tables

- * New Heavy Duty Actuators & Controller
- * New Programmable Hand Control (3 settings)
- * New Manual Foot Control (hands free operation)
- * Old wiring & hydraulic system completely removed
- * Warranty: Two (2) Years

\$3650.00
Installed on Site

No More Oil !!!

Revolutionary
NEW
Product



ELECTRAFIT
Advanced Technology Applied

stryker® Stryker 1001 Stretcher -

Refurbished

- 83" x 26.5"
- Back Up/Down (Pneumatic Assist)
- 8" Casters
- Fold-Under Siderails
- 4 Wheel Locking
- Steering
- 500 Lbs Weight Capacity
- 4" Pad (new)

\$2150 plus freight
1-year parts warranty



BLICKMAN

\$1895.00 + Shipping

ADJUSTABLE HEIGHT INSTRUMENT TABLE

Designed to straddle the Procedure table.
18 gauge, 304 stainless construction

Adjustment Height 36" to 56"
36" long
20" wide
3" Casters (two with brakes)

Cordell Medical

Call Today!

888-655-3534

Richmond, Virginia

www.CordellMed.com