

[the 950 series]



951



951
[black composite base]

Physician Seating Specifications

950 Series Exam Stool:

Back Rest Option



- single lever hand adjustment
- standard height adj.: 17" to 22"
- seat diameter: 15.75"
- base diameter: 24"
- high density fire resistant foam
- 3/4" thick composite substrate with steel reinforcement plate

Series Options

Back Rest: (specify "951-B")

- adjustable back
- seat back dimensions:
15.75" x 9.5"

Foot Rest: (specify "951-FR")

- polished chrome finish
- foot rest diameter: 18"

Height Adjustment:

- 19" to 26" - specify "951-7"
- 23" to 33" - specify "951-10"

18" Foot Rest Option




intensa

medical & laboratory seating

P.O. Box 5981 High Point, NC 27262
Toll-Free: 1-888-636-2660 Fax: (336) 884-4007

www.intensa.net