

BrewerAccess

EXAM TABLE



Advanced Safety

- Large patient step
- Stable four leg step design
- Largest step surface available: 17.7" x 20.4" (45 cm x 51.8 cm)

Patient Comfort

- Simple pneumatic backrest for easy patient positioning
- Multi-density, bolstered support upholstery
- Leg extension/plastic debris tray
- Adjustable ergonomic stirrups

Abundant Storage

- Patented pass-through side drawers accessible from right or left side of the table
- Drawer stops (provided) allow drawer to open from only one side if desired
- Front drawer: 19.8" x 19.2" x 3.8" (50.3 cm x 48.8 cm x 9.7 cm)
- Side drawers: 22.7" x 20.2" x 7" (57.4 cm x 51.3 cm x 17.8 cm)
- Additional paper storage under backrest for up to 6 additional paper rolls
- Paper roll storage capacity: 23" x 23" x 3.25" (58.4 cm x 58.4 cm x 8.3 cm)
- Slide-through paper roll holder accepts 21" & 18" rolls (53.3 cm & 45.7 cm)

Table Specifications

- 500 lb (226 kg) weight capacity
- Standard 5-year warranty
- Height: 32" (81.3 cm)
- Width: 28" (71.1 cm)
- Length: 57.3" extends to 71" (145.5 cm – 180.3 cm)
- Leg extension length: 15" (38.1 cm)
- Step Height: 8" (20.3 cm)
- Step Length: 17.7" (45 cm)
- Step Width: 20.4" (51.8 cm)
- Hospital-grade duplex electrical outlets
- 115VAC, 60Hz, 5.8 amps max
- Outlets can be put on left or right side
- Easy to clean seamless upholstery
- 12 standard upholstery color options
- Custom colors available
- Optional: 5° pelvic tilt/60-Watt drawer heater in front top drawer