



Assist Power Procedure Table

BrewerAssist POWER PROCEDURE TABLE

The Brewer **Assist Power Procedure Table™** can be lowered down to 19" making wheel chair transfers effortless. The 450 lb weight capacity and the additional 5" leg extension accommodates larger patients more comfortably for safer clinical outcomes. As with all Brewer Power Tables, the **Brewer Assist** comes with our standard 3 year warranty and 12 standard seamless upholstery colors.

Product Benefits

- 450 lb weight capacity
- Standard 3 year warranty
- Most technologically advanced
- Lowers to 19" for lateral patient transfer
- Exclusive three point pivot headrest
- Seamless easy to clean upholstery
- Footrest extends an additional 5"
- Intuitive foot controls
- **12 standard colors** with optional upgrades


8411 Glazebrook Avenue
Richmond VA 23228
804-262-0662