

Preventive Maintenance Checklist

COLPOSCOPE

PM Date: _____	Make: _____
Customer: _____	Model: _____
Facility: _____	Serial #: _____
Room: _____	Asset #: _____
	Tech: _____

A. CLEANING & ADJUSTMENTS

_____ 1. Oculars
 _____ Clean External
 _____ Clean Internal

_____ 2. Body
 _____ Clean External
 _____ Clean lens

_____ 3. Head Ways
 _____ External lens
 _____ Adjust

_____ 4. Ocular Tubes
 _____ Lube Adjustable
 _____ Clean External

_____ 5. Objectives
 _____ Clean External Lens
 _____ Parafoal

_____ 6. Coarse Adjustment
 _____ Lube Ways/Bearings
 _____ Adjust Tension
 _____ Lube Pinions

_____ 7. Light Source
 _____ Clean External
 _____ Center & Align
 _____ Clean Transformer
 _____ Check Controls

B. FINAL CHECK

_____ 1. Check Objectives w/ Prepared Slides
 _____ Clarity
 _____ Set Up
 _____ 2. Inspect Line Cord for Electrical Safety
 _____ 3. Customer Available to Check

C. OVERALL CONDITION

_____ Good
 _____ Fair
 _____ Poor

Technicians Notes: _____

This inspection was performed in accordance with manufacturer specifications.