

CERTIFIED PREVENTIVE MAINTENANCE PROCEDURE CHECKLIST

Urological Flow Meters

PM Date: _____	Make: _____
Customer: _____	Model: _____
Facility: _____	Serial #: _____
Room: _____	Asset #: _____
	Tech: _____

Description of weights:

10mg to 100g Polished Weight Set, ASTM Class 1, SN _____

Calibration Checklist:

_____	Linearity
_____	Repeatability
_____	Readability
_____	Weighing Capacity
_____	Taring Range
_____	Sensitivity Drift
_____	Stabilization Time

Weight Used:

_____	50 g
_____	100 g
_____	500 g
_____	1000 g
_____	2000 g

Weight Recorded:

_____	g	tolerance +/- 3
_____	g	tolerance +/- 5
_____	g	tolerance +/- 8
_____	g	tolerance +/- 15
_____	g	tolerance +/- 20

Operational Overview:

_____	Cleaned
_____	Adjusted

Overall Condition:

_____	Pass
_____	Fail/ Require Service